

Sharps Rules - Reuse for Multidosing

NEVER Recap Needles. EVER.

When reusing needles for multidosing:

Use an appropriate safety glide needle



Use safety glide mechanism when administration of doses are completed to cover needle. The cover cannot be removed after it is engaged.

Place safety glide needle in conical tube between doses



Use a conical tube on its side and place the syringe inside when not in use.

Sharps Rules - Replacing the Needle

NEVER Recap Needles. EVER.

Steps for Removing/Replacing Needles:

1. Use hemostats to grasp base of needle and twist off.

Required for ALL needles, even safety needles.

2. Place **IMMEDIATELY** into sharps container.
3. New, capped needle may be attached using hands.

